

## Potassium hydrogen sulphate

**GRM1790**

### Product Identifier

CAS No.	:	7646-93-7
EC No.	:	231-594-1
Molecular Formula	:	KHSO <sub>4</sub>
Molecular Weight	:	136.17
Synonym	:	Potassium bisulphate
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Assay	:	97.00 - 102.00%

### Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H314- H335
Precautionary Statement(s)	:	P261- P280- P305+P351+P338- P310
UN No.	:	2509
Class	:	8
Packing Group	:	II
RTECS	:	TS7200000
WGK	:	1