

Dithizone, Hi-AR®/ACS

GRM1816


Product Identifier

CAS No.	:	60-10-6
EC No.	:	200-454-1
Molecular Formula	:	C ₁₃ H ₁₂ N ₄ S
Molecular Weight	:	256.33
HS Code	:	2930 90 99
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Purple to purple-black to black crystals or powder or chunks or solid
Solubility	:	Insoluble in water
Heavy metals (as Pb)	:	<= 0.002%
Residue on ignition	:	<= 0.30%
Ratio of absorbances	:	>= 1.55
Assay (UV)	:	min. 85.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319
Precautionary Statement(s)	:	P280- P302 + P352- P305+P351+P338- P332 + P313- P337 + P313- P362 + P364
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LQ9450000
WGK	:	1