

## Eschka's mixture, Hi-AR<sup>®</sup>

**GRM1819**

(For determination of Sulphur in coal)


### Product Identifier

CAS No.	:	8007-09-8
Molecular Formula	:	2MgO.Na <sub>2</sub> CO <sub>3</sub>
HS Code	:	3822 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White powder
Solubility	:	33.3 mg soluble in 1 mL of dilute hydrochloric acid
Hydrochloric acid insoluble matter	:	<= 0.005%
Sulphate & Sulphite (as SO <sub>4</sub> )	:	<= 0.005%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P280- P305+P351+P338- P337 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1