

GRM1846 Salol**Product Number Packing**

GRM1846	: 100G
GRM1846	: 500G


Product Information

Product Code	: GRM1846
Product Name	: Salol
Synonym	: Phenyl salicylate
Molecular Formula	: C ₁₃ H ₁₀ O ₃
Molecular Weight	: 214.22
CAS No.	: 118-55-8
EC No.	: 204-259-2
HS Code	: 2918 23 00
Shelf Life	: 4 years


Technical Specification

Appearance	: White crystals or powder or granules or chunks
Solubility	: 33.3 mg soluble in 1 mL of methanol
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 41 - 43°C
Assay (GC)	: min. 98.00 %

GHS Safety Information

Hazard Statement(s)	: H315-H319-H335
Precautionary Statement(s)	: P261-P305+P351+P338
Signal Word	: Warning
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 36/37/38
S-Phrase(s)	: 26-36
WGK	: 2
RTECS	: VO6125000
Flash Point(°F)	: 235.4 °F
Flash Point(°C)	: 113 °C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

Transport Information

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein



Product Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein