

## Sodium (meta)arsenite, Hi-AR™

**GRM1847**


### Product Identifier

CAS No.	:	7784-46-5
EC No.	:	232-070-5
EC Index No	:	033-005-00-1
Molecular Formula	:	NaAsO <sub>2</sub>
Molecular Weight	:	129.91
HS Code	:	2842 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to yellow to grey powder.
Solubility	:	33.3 mg soluble in 1 mL of water to give clear to turbid solution
Assay (Iodometry)	:	99.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300+H310- H331- H350- H410
Precautionary Statement(s)	:	P201- P261- P264- P273- P280- P301 + P310
UN No.	:	2027
Class	:	6.1
Packing Group	:	II
RTECS	:	CG3675000
WGK	:	3