
Acetanilide, Hi-AR™

GRM1865


Product Identifier

CAS No.	:	103-84-4
EC No.	:	203-150-7
Molecular Formula	:	C ₈ H ₉ NO
Molecular Weight	:	135.16
Synonym	:	N-Phenylacetamide
HS Code	:	2924 29 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown or grey crystals or powder or chips or flakes or chunks
Solubility	:	100 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	111 - 115°C
Assay (GC)	:	min. 99.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302
Precautionary Statement(s)	:	P301+P312+P330
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AD7350000
WGK	:	1