

Methylamine hydrochloride

GRM1881


Product Identifier

CAS No.	:	593-51-1
EC No.	:	209-795-0
Molecular Formula	:	CH ₅ N.HCl
Molecular Weight	:	67.52
Synonym	:	Methylammonium chloride
HS Code	:	2921 11 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to light yellow crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	226 - 236°C
Assay (AT)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	PA0603000
WGK	:	3