

## Methyl paraben

GRM1899


### Product Identifier

CAS No.	:	99-76-3
EC No.	:	202-785-7
Molecular Formula	:	C <sub>8</sub> H <sub>8</sub> O <sub>3</sub>
Molecular Weight	:	152.15
Synonym	:	Methyl-4-hydroxybenzoate
HS Code	:	2918 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Colourless to off-white crystals or powder
Solubility	:	100 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	125 - 128°C
Chloride (Cl)	:	<= 0.005%
Sulfate (SO <sub>4</sub> )	:	No turbidity should be produced within 10 minutes
Assay (NaOH Titration/HPLC)	:	99.00 - 101.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH2450000
WGK	:	1