

Antimony (III) sulphate

GRM1946


Product Identifier

CAS No.	:	7446-32-4
EC No.	:	231-207-6
Molecular Formula	:	Sb ₂ (SO ₄) ₃
Molecular Weight	:	531.71
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder or lumps
Solubility	:	Insoluble in water
Assay (Iodimetry/ICP)	:	97.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302+H332- H411
Precautionary Statement(s)	:	P261- P264- P270- P273- P301 + P312- P304 + P340
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	CC8225000
WGK	:	2