

---

## Bees wax for G.C.

**GRM1959**

---

### Product Identifier

CAS No.	:	8012-89-3
EC No.	:	232-383-7
HS Code	:	1521 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to pale yellow to yellow-tan solid or pieces or plates.
Solubility	:	33.3 mg soluble in 1 mL of chloroform (Clear to hazy)
Saponification value	:	87 - 104

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1