
Starch, Corn

GRM198

Product Identifier

CAS No.	:	9005-25-8
EC No.	:	232-679-6
Molecular Formula	:	(C ₆ H ₁₀ O ₅) _n
Molecular Weight	:	(162.14) _n
Synonym	:	Maize starch
HS Code	:	1108 12 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to slightly yellowish powder or irregular white masses which are readily reducible to powder
Solubility	:	Insoluble in water
Loss on drying (at 105°C, 1 hr)	:	<= 15.00%
Sulphated ash	:	<= 0.60%
Iron (Fe)	:	<= 0.004%
Oxidising substances	:	Passes test

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	1