

p-Bromophenol

GRM2012


Product Identifier

CAS No.	:	106-41-2
EC No.	:	203-394-4
Molecular Formula	:	C ₆ H ₅ BrO
Molecular Weight	:	173.01
Synonym	:	4-Bromophenol
HS Code	:	2908 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to orange to brown crystals or powder or chunks.
Solubility	:	100 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern.
Assay (GC)	:	>= 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315
Precautionary Statement(s)	:	P270- P280- P301 + P312- P302 + P352- P332 + P313- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SJ7960000
WGK	:	1