
Carbopol 940

GRM2033

Product Identifier

HS Code	:	3906 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White hygroscopic powder
Solubility	:	1% w/v aqueous solution gives a slight opalescent gel at room temperature.
Loss on drying (at 80°C, 2 hr)	:	<= 3.00%
Viscosity	:	40000 - 60000 cPs (0.5% neutral solution, spindle no. 7, 20 rpm at 25°C)
Assay (Carboxylic acid, NaOH Titration)	:	56.00 - 68.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1