



**GRM206**

**Soyabean meal, Defatted**

---

**Product Number**

**Packing**

GRM206 : 500G

---

**Product Information**

Product Code : GRM206  
Product Name : Soyabean meal, Defatted  
Shelf Life : 4 years

---

**Technical Specification**

Appearance : Cream coloured free flowing powder having characteristic odour.  
Solubility : Insoluble in water  
Loss on drying (at 105°C, 1 hr) : <= 8.00%  
Fat content : <= 1.00%  
Total nitrogen (On dry basis) : min. 8.00%  
Total protein (on dry basis) : min. 50.00%

---

**Risk and Safety Information**

WGK : 3  
Storage Temperature(°C) : Store below 30°C

---

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods