

**GRM2189 Ethylparaben****Product Number Packing**

GRM2189 : 250G

**Product Information**

Product Code : GRM2189  
Product Name : Ethylparaben  
Synonym : Ethyl-4-hydroxybenzoate; p-Hydroxybenzoic acid ethylester  
Molecular Formula : C<sub>9</sub>H<sub>10</sub>O<sub>3</sub>  
Molecular Weight : 166.17  
CAS No. : 120-47-8  
EC No. : 204-399-4  
HS Code : 2918 29 00  
Shelf Life : 4 years

**Technical Specification**

Appearance : White or almost white crystals or powder  
Solubility : 33.3 mg soluble in 1 mL of alcohol  
FTIR (KBr disc) : Matches with the standard pattern  
Melting range : 114 - 118°C  
Sulphated ash : <= 0.1%  
Assay (GC/NaOH Titration) : 99.00 - 102.00%

**GHS Safety Information**

Hazard Statement(s) : H315-H319-H335  
Precautionary Statement(s) : P261-P305+P351+P338  
Signal Word : Warning  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 36/37/38  
S-Phrase(s) : 26-36  
WGK : 1  
RTECS : DH2190000  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods



## Product Information

IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

---

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein