
(Hydroxypropyl) methyl cellulose

GRM2257

Product Identifier

CAS No.	:	9004-65-3
HS Code	:	3912 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellowish white fibrous or granular powder, hygroscopic afterdrying
Solubility	:	33.3 mg soluble in 1 mL of cold water
pH (1% in water at 25°C)	:	5.50 - 8.00
Identification test with tannic acid	:	Passes test
Loss on drying (at 105°C, 2 hr)	:	<= 10.00%
Sulphated ash	:	<= 3.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1