

Celite®

GRM226


Product Identifier

Synonym : Filter aid powder
HS Code : 2512 00 90
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : White to tan to grey powder
Solubility : Insoluble in water
Composition : SiO₂: 90.20%
Loss on ignition: 0.20%
Al₂O₃: 3.80%
Fe₂O₃: 1.20%
P₂O₅: 0.20%
TiO₂: 0.20%
CaO: 0.50%
MgO: 0.60%
Na₂O+K₂O: 3.10%

Safety Information

Hazard Pictogram(s) : 

Signal Word : Warning
Hazard Statement(s) : H319- H335- H373
Precautionary Statement(s) : P260- P280- P304 + P340- P305+P351+P338- P403 + P233- P501
UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1