

4,4'-Isopropylidenediphenol

GRM2273


Product Identifier

CAS No.	:	80-05-7
EC No.	:	201-245-8
EC Index No	:	604-030-00-0
Molecular Formula	:	C ₁₅ H ₁₆ O ₂
Molecular Weight	:	228.29
Synonym	:	Bisphenol A; 2,2-Bis(4-Hydroxyphenyl)propane
HS Code	:	2907 23 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White powder or crystals or granules or flakes or chips or beads
Solubility	:	33.3 mg soluble in 1 mL of methanol
FTIR	:	Matches with the standard pattern
Melting range	:	150 - 160°C
Assay (GC)	:	>= 97.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H317- H318- H335- H360F- H410
Precautionary Statement(s)	:	P261- P273- P280- P305+P351+P338- P362 + P364- P403 + P233
EU Phrases	:	EUH380- EUH430
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	SL6300000
WGK	:	3