
Manganese (II) acetate tetrahydrate, Hi-AR™

GRM2298


Product Identifier

CAS No.	:	6156-78-1
EC No.	:	211-334-3
Molecular Formula	:	C ₄ H ₆ MnO ₄ .4H ₂ O
Molecular Weight	:	245.09
Synonym	:	Manganous acetate
Storage	:	Below 30°C
Shelf life	:	4 Years

Technical Specification

Appearance	:	White to faint red crystals
Solubility	:	100 mg soluble in 1 mL of water
Chloride (Cl)	:	<= 0.001%
Heavy metals (as Pb)	:	<= 0.0005%
Iron (Fe)	:	<= 0.0005%
Sulfate (SO ₄)	:	<= 0.005%
Assay (EDTA Titration)	:	99.50 - 101.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AI5775000
WGK	:	2