

---

## Fuller's earth

**GRM232**

---

### Product Identifier

CAS No.	:	8031-18-3
HS Code	:	2508 40 90
Storage	:	Below 30°C
Shelf life	:	4 years


---

### Technical Specification

Appearance	:	Pale yellow to brown powder
Solubility	:	Insoluble in water
Particle size 1	:	min. 80% passing through 150 mesh (106 µm)
Particle size 2	:	max. 20% retain on 250 mesh (63 µm)

---

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H332- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	3