

---

## Methyl salicylate

**GRM2331**

---

### Product Identifier

CAS No.	:	119-36-8
EC No.	:	204-317-7
Molecular Formula	:	C <sub>8</sub> H <sub>8</sub> O <sub>3</sub>
Molecular Weight	:	152.15
Synonym	:	Methyl-2-hydroxybenzoate
HS Code	:	2918 23 10
Storage	:	Below 30°C
Shelf life	:	4 years


---

### Technical Specification

Appearance	:	Colourless to light yellow liquid
Solubility	:	1 mL miscible in 1 mL of ethanol
FTIR	:	Matches with the standard pattern
Refractive index (n <sub>20/D</sub> )	:	1.534 - 1.538
Density (at 25°C)	:	1.175 - 1.185 g/mL
Assay (GC/NaOH Titration)	:	99.00 - 100.50%

---

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VO4725000
WGK	:	1