

**GRM2363 Palladium (II) acetate****Product Number Packing**

GRM2363 : 1G

**Product Information**

Product Code : GRM2363  
Product Name : Palladium (II) acetate  
Synonym : Palladium diacetate; Acetic acid palladium (II) salt  
Molecular Formula : C<sub>4</sub>H<sub>6</sub>O<sub>4</sub>P d  
Molecular Weight : 224.51  
CAS No. : 3375-31-3  
EC No. : 222-164-4  
HS Code : 2843 90 90  
Shelf Life : 4 years

**Technical Specification**

Appearance : Light orange to brown powder  
Solubility : 33.3 mg soluble in 1 mL of chloroform  
Assay (Gravimetry) : 99.00 - 102.50%

**GHS Safety Information**

Hazard Statement(s) : H318  
Precautionary Statement(s) : P280-P305+P351+P338  
Signal Word : Danger  
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 41  
S-Phrase(s) : 26-39  
WGK : 2  
RTECS : AJ1900000  
Storage Temperature(°C) : Store below 30°C  
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods