
Pectinase, From *Aspergillus niger*

GRM2366

Product Identifier

| | | |
|------------|---|------------|
| CAS No. | : | 9032-75-1 |
| EC No. | : | 232-885-6 |
| HS Code | : | 3507 90 99 |
| Storage | : | 2 - 8°C |
| Shelf life | : | 3 years |

Technical Specification

| | | |
|---------------------------------|---|--|
| Appearance | : | White to pale amber powder or solid |
| Solubility | : | 10 mg soluble in 10 mL water |
| Loss on drying (at 105°C, 1 hr) | : | <= 10.00% |
| Enzyme activity | : | 8000 - 12000 U/g |
| Unit definition | : | One unit will liberate 1.0 micromole of galacturonic acid from polygalact-uronic acid per minute at pH 4.0 at 25°C |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 3 |