
Phthalimide potassium

GRM2387


Product Identifier

CAS No.	:	1074-82-4
EC No.	:	214-046-6
Molecular Formula	:	C ₈ H ₄ KNO ₂
Molecular Weight	:	185.22
Synonym	:	Potassium phthalimide
HS Code	:	2925 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow to green crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Assay (NT)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	3