

## Selenium dioxide

**GRM2434**


### Product Identifier

CAS No.	:	7446-08-4
EC No.	:	231-194-7
EC Index No	:	034-002-00-8
Molecular Formula	:	SeO <sub>2</sub>
Molecular Weight	:	110.96
HS Code	:	2811 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to pink to beige crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
Assay (ICP/Gravimetry)	:	min. 98.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301+H331- H373- H410
Precautionary Statement(s)	:	P260- P273- P280- P301 + P310- P361 + P364- P501
UN No.	:	3283
Class	:	6.1
Packing Group	:	III
RTECS	:	VS8575000
WGK	:	2