

Selenium sulfide

GRM2435


Product Identifier

CAS No.	:	7488-56-4
EC No.	:	231-303-8
Molecular Formula	:	SeS ₂
Molecular Weight	:	143.09
HS Code	:	2813 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Bright orange to orange to reddish-brown powder
Solubility	:	Insoluble in water
Assay (% Se, Na ₂ S ₂ O ₃ T/Gravimetry)	:	52.00 - 55.50 %

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301- H331- H373- H410
Precautionary Statement(s)	:	P261- P273- P301 + P310- P311- P501
UN No.	:	2657
Class	:	6.1
Packing Group	:	II
RTECS	:	VS8925000
WGK	:	3