
Leishman's stain, Hi-Cert®

GRM244

Product Identifier

CAS No.	:	12627-53-1
EC No.	:	235-732-1
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Dark green to black powder
Solubility	:	10 ppm in methanol yields blue solution with greenish fluorescence
Suitability for microscopy	:	Passes test
Absorption maxima 1	:	640 - 650 (10 ppm in methanol)
Absorption maxima 2	:	515 - 530 (10 ppm in methanol)
Absorptivity 1	:	>= 650 (10 ppm in methanol at 513 - 530 nm)
Absorptivity 2	:	>= 900 (10 ppm in methanol at 640 - 650 nm)
Loss on drying (at 105°C, 1 hr)	:	<= 8.0%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1