

Sodium tetrafluoroborate

GRM2445


Product Identifier

CAS No.	:	13755-29-8
EC No.	:	237-340-6
Molecular Formula	:	BF ₄ Na
Molecular Weight	:	109.79
Synonym	:	Sodium fluoroborate
HS Code	:	2826 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown crystals or powder or needles or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Titration)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P501
UN No.	:	3260
Class	:	8
Packing Group	:	II
RTECS	:	ED2975000
WGK	:	1