

GRM2448 Tin (IV) chloride pentahydrate (Stannic chloride pentahydrate)**Product Number Packing**

GRM2448 : 500G

Product Information

Product Code : GRM2448
Product Name : Tin (IV) chloride pentahydrate
Synonym : Stannic chloride pentahydrate
Molecular Formula : $\text{SnCl}_4 \cdot 5\text{H}_2\text{O}$
Molecular Weight : 350.60
CAS No. : 10026-06-9
HS Code : 2827 39 90
Shelf Life : 4 years

Technical Specification


Appearance : White to yellow crystals or powder or chunks or lumps
Solubility : 33.3 mg soluble in 1 mL of water
Iron (Fe) : $\leq 0.005\%$
Sulfate (SO_4) : $\leq 0.02\%$
Assay (AT/ICP) : 98.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H314
Precautionary Statement(s) : P280-P305+P351+P338-P310
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 34
S-Phrase(s) : 26-36/37/39-45
WGK : 3
RTECS : XP8870000
Storage Temperature($^{\circ}\text{C}$) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 2440
Class : 8
Packaging Group : 3
Marine Pollutant : No



Product Information

ADR/RID : 2440 8/PG 3
IMDG : 2440 8/PG 3
IATA : 2440 8/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein