

---

## Tetrabutylammonium hydroxide, 40% solution in water

**GRM2463**

### Product Identifier

---

CAS No.	:	2052-49-5
EC No.	:	218-147-6
Molecular Formula	:	C <sub>16</sub> H <sub>37</sub> NO
Molecular Weight	:	259.47
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years


### Technical Specification

---

Appearance	:	Colourless to yellow liquid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (HCl Titration)	:	38.00 - 42.00%

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P310
UN No.	:	3267
Class	:	8
Packing Group	:	II
WGK	:	1