

Tetrabutylammonium hydroxide 40% solution in methanol

GRM2464


Product Identifier

CAS No.	:	2052-49-5
EC No.	:	218-147-6
Molecular Formula	:	C ₁₆ H ₃₇ NO
Molecular Weight	:	259.47
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colorless to yellow liquid
Solubility	:	1 mL miscible in 1 mL of methanol
Assay (HCl Titration)	:	37.0 - 43.0%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H225- H301+H311+H331- H314- H370
Precautionary Statement(s)	:	P210- P260- P280- P301 + P310- P305+P351+P338- P303+P361+P353
UN No.	:	3286
Class	:	3 (6.1,8)
Packing Group	:	II
RTECS	:	Not available
WGK	:	2