

---

## Nigrosine

**GRM247**

---

### Product Identifier

|            |   |            |
|------------|---|------------|
| CAS No.    | : | 8005-03-6  |
| C.I. No.   | : | 50420      |
| HS Code    | : | 3204 19 90 |
| Storage    | : | Below 30°C |
| Shelf life | : | 4 years    |

---

### Technical Specification

|                                 |   |  |
|---------------------------------|---|--|
| Appearance                      | : | Brown to black crystals or powder          |
| Solubility                      | : | 25 ppm in 50% ethanol gives clear solution |
| Absorption maxima               | : | 570 - 580 nm (25 ppm in 50% ethanol)       |
| Loss on drying (at 110°C, 2 hr) | : | <= 15.00%                                  |

---

### Safety Information

|               |   |                     |
|---------------|---|---------------------|
| UN No.        | : | Not dangerous goods |
| Class         | : | -                   |
| Packing Group | : | -                   |
| RTECS         | : | GC4762250           |
| WGK           | : | 1                   |