

GRM2486 p-Toluenesulphonyl chloride**Product Number Packing**

GRM2486 : 500G

Product Information

Product Code : GRM2486
Product Name : p-Toluenesulphonyl chloride
Synonym : Tosyl chloride
Molecular Formula : C₇H₇ClO₂S
Molecular Weight : 190.65
CAS No. : 98-59-9
EC No. : 202-684-8
HS Code : 2904 90 95
Shelf Life : 4 years

Technical Specification


Appearance : White to yellow solid or crystals or powder or chunks or granules
Solubility : 100 mg soluble in 1 mL of alcohol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 65 - 69°C
Assay (AT/GC) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H315-H318
Precautionary Statement(s) : P280-P305+P351+P338
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 38-41
S-Phrase(s) : 26-39
WGK : 1
RTECS : DB8929000
Flash Point(°F) : 262.4°F
Flash Point(°C) : 128°C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 3261



Product Information

Class	:	8
Packaging Group	:	3
Marine Pollutant	:	No
ADR/RID	:	3261 8/PG 3
IMDG	:	3261 8/PG 3
IATA	:	3261 8/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein