

GRM2487 o-Toluic acid

Product Number Packing

GRM2487 : 500G

Product Information

Product Code : GRM2487
Product Name : o-Toluic acid
Synonym : 2-Methylbenzoic acid
Molecular Formula : C₈H₈O₂
Molecular Weight : 136.15
CAS No. : 118-90-1
EC No. : 204-284-9
HS Code : 2916 39 90
Shelf Life : 4 years

Technical Specification

Appearance : White to yellow-tan needle like crystals or flakes or powder
Solubility : 33.3 mg soluble in 1 mL of alcohol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 101 - 107°C
Assay (NaOH Titration/GC) : 98.00 - 102.00%


GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36
WGK : 3
RTECS : XU1400000
Flash Point(°F) : 298.4°F
Flash Point(°C) : 148 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No



Product Information

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein