

GRM2512 Urethane**Product Number Packing**

GRM2512	: 100G
GRM2512	: 500G

Product Information

Product Code	: GRM2512
Product Name	: Urethane
Synonym	: Ethylcarbamate
Molecular Formula	: C ₃ H ₇ O ₂ N
Molecular Weight	: 89.09
CAS No.	: 51-79-6
EC No.	: 200-123-1
HS Code	: 2924 19 00
Shelf Life	: 4 years

Technical Specification

Appearance	: Colourless crystals or chunks
Solubility	: 33.3 mg soluble in 1 mL of water
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 45 - 55°C
Assay (Total nitrogen/GC)	: 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s)	: H302-H350
Precautionary Statement(s)	: P201-P308+P313
Signal Word	: Danger
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 45
S-Phrase(s)	: 53-45
WGK	: 3
RTECS	: FA8400000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

**Transport Information**

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods



Product Information

IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods