

Wright's stain, Hi-Cert[™]

GRM265


Product Identifier

CAS No.	:	68988-92-1
EC No.	:	273-541-5
HS Code	:	3204 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Dark green to purple to very dark blue to black solid or powder
Solubility	:	10 ppm in methanol yields clear solution
Absorption maxima 1	:	640 - 648 nm (10 ppm in methanol)
Absorption maxima 2	:	522 - 525 nm (10 ppm in methanol)
Loss on drying (at 110°C, 2 hr)	:	<= 10.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H318
Precautionary Statement(s)	:	P280- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1