

GRM2670 4'-Bromoacetophenone**Product Number Packing**

GRM2670	: 100G
GRM2670	: 1KG
GRM2670	: 250G

Product Information

Product Code	: GRM2670
Product Name	: 4'-Bromoacetophenone
Synonym	: p-Bromoacetophenone; 1-Acetyl-4-bromobenzene
Molecular Formula	: C ₈ H ₇ BrO
Molecular Weight	: 199.04
CAS No.	: 99-90-1
EC No.	: 202-799-3
HS Code	: 2914 70 00
Shelf Life	: 4 years

Technical Specification

Appearance	: White to yellow crystals or flakes or powder
Solubility	: 33.3 mg soluble in 1 mL of alcohol
FTIR (KBr disc)	: Matches with the standard pattern
Melting range	: 48 - 54°C
Assay (GC)	: min. 98.00%

GHS Safety Information

Hazard Statement(s)	: H315-H319-H334-H335
Precautionary Statement(s)	: P261-P305+P351+P338-P342+P311
Signal Word	: Danger
Hazard Pictogram(s)	

**Risk and Safety Information**

R-Phrase(s)	: 22-36/37/38-42/43
S-Phrase(s)	: 22-26-36/37
WGK	: 3
Flash Point(°F)	: 235.4°F
Flash Point(°C)	: 113°C
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	





Product Information

Transport Information

Marine Pollutant	:	No
ADR/RID	:	Not Dangerous Goods
IMDG	:	Not Dangerous Goods
IATA	:	Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein