

Carbazole

GRM2696


Product Identifier

CAS No.	:	86-74-8
EC No.	:	201-696-0
Molecular Formula	:	C ₁₂ H ₉ N
Molecular Weight	:	167.21
HS Code	:	2933 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown crystals or powder or chunks or flakes
Solubility	:	33.3 mg soluble in 1 mL of acetone
FTIR	:	Matches with the standard pattern
Melting range	:	240 - 247°C
Assay (GC)	:	min. 95.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H319- H351- H373- H410
Precautionary Statement(s)	:	P201- P260- P273- P280- P308 + P313- P314
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	FE3150000
WGK	:	2