
5-Methyluracil

GRM276

Product Identifier

CAS No.	:	65-71-4
EC No.	:	200-616-1
Molecular Formula	:	C ₅ H ₆ N ₂ O ₂
Molecular Weight	:	126.11
Synonym	:	Thymine
HS Code	:	2933 59 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white crystals or powder or solid
Solubility	:	10 mg soluble in 10 mL of DMSO
FTIR	:	Matches with the standard pattern
Assay (HPLC)	:	min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XP2100000
WGK	:	2