

N,N-Dimethylformamide dimethyl acetal

GRM2788


Product Identifier

CAS No.	:	4637-24-5
EC No.	:	225-063-3
Molecular Formula	:	C ₅ H ₁₃ NO ₂
Molecular Weight	:	119.16
HS Code	:	2924 19 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to light yellow liquid
Solubility	:	1 mL miscible in 1 mL of water
FTIR	:	Matches with the standard pattern
Refractive index (n 20/D)	:	1.391 - 1.401
Density (at 25°C)	:	0.892 - 0.902 g/mL
Assay (GC/NT)	:	95.00 - 107.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H225- H315- H317- H319- H332- H335
Precautionary Statement(s)	:	P210- P261- P280- P304+P340+P312- P305+P351+P338- P370 + P378
UN No.	:	3271
Class	:	3
Packing Group	:	II
RTECS	:	Not available
WGK	:	1