
Diphenylcarbazone, Hi-AR™/ACS (Mixture with s-diphenylcarbazine)

GRM2797

Product Identifier

CAS No.	:	538-62-5 with CAS No. : 140-22-7
EC No.	:	208-698-0
Molecular Formula	:	C ₁₃ H ₁₂ N ₄ O
Molecular Weight	:	240.26
HS Code	:	2928 00 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Orange to brownish orange crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of acetone
Sensitivity	:	Passes
Melting range	:	154 - 161°C
Residue after ignition	:	<= 0.1%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	LQ9420000
WGK	:	3