
Benzethonium chloride (0.004M solution)

GRM2845

Product Identifier

CAS No.	:	121-54-0
Molecular Formula	:	C ₂₇ H ₄₂ ClNO ₂
Molecular Weight	:	448.08
Synonym	:	Hyamine [®] 1622
HS Code	:	2923 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless solution
Solubility	:	33.3 mg soluble in 1 mL of water
Concentration	:	0.003996 - 0.004004 mol/L

Safety Information

Hazard Statement(s)	:	H412
Precautionary Statement(s)	:	P273- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1