

Hydroquinone dimethyl ether

GRM2849


Product Identifier

CAS No.	:	150-78-7
EC No.	:	205-771-9
Molecular Formula	:	C ₈ H ₁₀ O ₂
Molecular Weight	:	138.16
HS Code	:	2909 30 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to faint beige deliquescent crystals or powder or chunks or chips or solid
Solubility	:	33.3 mg soluble in 1 mL methanol
FTIR	:	Matches with the standard pattern
Melting range	:	54 - 59°C
Assay (GC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CZ6650000
WGK	:	1