
Kieselguhr, Hi-AR[®]

GRM2887


Product Identifier

CAS No.	:	91053-39-3
HS Code	:	2512 00 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white powder
Solubility	:	Insoluble in water
pH (10% suspension)	:	5.00 - 10.00
Acid soluble Iron (Fe)	:	<= 0.04%
Acid soluble Sulfate (SO ₄)	:	<= 0.02%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H335
Precautionary Statement(s)	:	P261- P271- P280- P304 + P340- P305+P351+P338- P312
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1