

Ampyrone, Hi-ARTM/ACS

GRM289


Product Identifier

CAS No.	:	83-07-8
EC No.	:	201-452-3
Molecular Formula	:	C ₁₁ H ₁₃ N ₃ O
Molecular Weight	:	203.24
Synonym	:	4-Aminoantipyrine, Free base, Hi-AR TM /ACS
HS Code	:	2933 11 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to yellow-brown crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	105 - 110°C
Assay (NT/HPLC)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CD2480000
WGK	:	1