

Phenacetin

GRM2941


Product Identifier

CAS No.	:	62-44-2
EC No.	:	200-533-0
Molecular Formula	:	C ₁₀ H ₁₃ NO ₂
Molecular Weight	:	179.22
HS Code	:	2924 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder or powder with lumps
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	133 - 138°C
Assay (GC/HPLC)	:	min. 98.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H350
Precautionary Statement(s)	:	P201- P202- P264.- P270- P301 + P312- P308 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AM4375000
WGK	:	3