

---

## Potassium pyroantimonate trihydrate, Hi-AR™

**GRM2952**

### Product Identifier

---

CAS No.	:	12208-13-8
EC No.	:	235-387-7
Molecular Formula	:	KSb(OH) <sub>6</sub>
Molecular Weight	:	262.90
HS Code	:	2841 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to off-white crystals powder
Solubility	:	10 mg soluble in 10 mL of boiling water gives hazy solution
Assay (Iodometry, on ignited basis)	:	99.00 - 102.00%

### Safety Information

---

Hazard Pictogram(s) :



Signal Word	:	Warning
Hazard Statement(s)	:	H302+H332- H411
Precautionary Statement(s)	:	P273
UN No.	:	1549
Class	:	6.1
Packing Group	:	III
WGK	:	2