

Antimony trioxide, Hi-LR™

GRM3038


Product Identifier

CAS No.	:	1309-64-4
EC No.	:	215-175-0
Molecular Formula	:	Sb ₂ O ₃
Molecular Weight	:	291.52
HS Code	:	2825 80 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder or solid
Solubility	:	10 mg soluble in 10 mL of dilute hydrochloric acid
Assay (Iodimetry)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H351
Precautionary Statement(s)	:	P201- P280- P308 + P313- P405- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	CC5650000
WGK	:	1