
Chromium (III) potassium sulphate dodecahydrate, Hi-LR™ GRM3043

Product Identifier

CAS No.	:	7788-99-0
EC No.	:	233-401-6
Molecular Formula	:	$\text{CrK}(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$
Molecular Weight	:	499.40
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Purple to violet to black crystals or powder or solid or chunks
Solubility	:	100 mg soluble in 1 mL of water
Chloride (Cl)	:	$\leq 0.03\%$
Iron (Fe)	:	$\leq 0.03\%$
Assay (Iodometry)	:	98.50 - 102.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GB6850000
WGK	:	2