

---

## Rubidium chloride, Hi-AR®

**GRM3060**

### Product Identifier

---

CAS No.	:	7791-11-9
EC No.	:	232-240-9
Molecular Formula	:	RbCl
Molecular Weight	:	120.92
HS Code	:	2827 39 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White crystals or powder or chunks or solid
Solubility	:	100 mg soluble in 1 mL of water
Magnesium (Mg)	:	<= 0.0005%
Heavy metals (as Pb)	:	<= 0.0003%
Iron (Fe)	:	<= 0.0003%
Sulphate (SO <sub>4</sub> )	:	<= 0.005%
Assay (AT back Titration)	:	min. 99.50%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VL8575000
WGK	:	1